

MEDUMAT Easy^{CPR}

The Guiding Emergency Ventilator







MEDUMAT Easy CPR



Safe and Effective Ventilation for Non-Intubated and Intubated Patients

Airway management and ventilation are some of the most important, yet challenging situations EMS faces in the pre-hospital setting. And ventilation continues to be one of the most difficult aspects of high-performance CPR. With MEDUMAT Easy^{CPR}, ventilation is consistent and reliable, lowering the risk of hyperventilation and gastric insufflation that comes with

manual bag-valve mask ventilation. MEDUMAT Easy^{CPR} can be used with a basic mask, a supraglottic airway or an endotracheal tube, allowing responders at all levels to maintain appropriate respiratory rates and tidal volumes throughout patient care.

The ventilator ensures the highest degree of safety by giving clear audio responses that guide the EMS responder through CPR procedures. When developing MEDUMAT Easy^{CPR}, we concentrated on ergonomic design, intuitive operation and a size and weight that would be easy to transport. Mounted on a LIFE-BASE portable system MEDUMAT Easy^{CPR} and our defibrillator/monitor MEDUCORE Standard² make the perfect emergency team.

Main Benefits

- Can be integrated easily in the standard process for CPR
- Simple to use with intuitive controls and predefined settings
- Small and lightweight
- Allows one person to administer two mechanical breaths in less than 5 sec.
- Can be used on one portable system together with MEDUCORE Standard or MEDUCORE Standard²
- Voice assistance available in 20 languages
- Reliable rate and volume not too much, not too little
- Decreased risk of hyperventilation and gastric insufflation compared to a BVM
- Demand flow mode for oxygenation of patients who are spontaneously breathing



Simple and Intuitive

MEDUMAT Easy^{CPR}'s built-in adult and pediatric settings allow you to quickly ventilate at the right rate and volume during critical situations. With minimal training, paramedics of all levels can use MEDUMAT Easy^{CPR} to better ventilate patients with a variety of conditions. And at less than 1 kg, MEDUMAT Easy^{CPR} is ideal for the pre-hospital setting.



Product Features

- 1. CPR mode for applying manual breaths via MEDUtrigger during resuscitation
- 2. One single knob for setting tidal volume and ventilation rate
- 3. Color coding for simple settings by patient age
- 4. Built-in visual and audible alarm system
- 5. Demand-flow mode supports spontaneous ventilation
- 6. Two pressure limits for tube and mask ventilation
- 7. Continuous monitoring of airway pressure via LED display
- 8. Connection for oxygen hose
- 9. Loudspeaker for voice assistance and metronome (can be deactivated)

More Consistent Mask Ventilation

Using the integrated MEDUtrigger, a single EMS professional can maintain a good mask seal with both hands while simultaneously ventilating the patient. Located on the ventilation mask, the MEDUtrigger allows the rescuer to deliver a single mechanical breath at a fixed tidal volume – simply by pressing a button with their thumb. Once a supraglottic airway or endotracheal tube is in place, MEDUMAT Easy^{CPR} can be used for continuous ventilation with a fixed ventilation rate and volume. This allows the EMS professional to accomplish other tasks besides ventilation.





炒 LIFE-BASE Portable Systems



Mobility Concept for MEDUMAT Easy^{CPR}

Available on LIFE-BASE portable systems, MEDUMAT Easy^{CPR} offers you many benefits in ventilation and oxygen therapy.

- 1. LIFE-BASE 3 NG with MEDUCORE Standard² and MEDUMAT Easy^{CPR} (WM 9932)
- Combination options with monitor/defibrillator
 MEDUCORE Standard² or
 MEDUCORE Standard (not shown)
- 3. Easy to transport
- 4. Both ventilation and oxygen inhalation possible in one system
- 5. Practical and safe storage of accessories



LIFE-BASE light with MEDUMAT Easy^{CPR} and MODUL CPAP (WM 28235)



LIFE-BASE III with MEDUMAT Easy^{CPR} (WM 9155)



LIFE-BASE III MASCAL
with MEDUMAT Easy^{CPR},
MODUL Oxygen and 3-way distributor
(WM 8000)

WEINMANN medical technology

10 Accessories and Replacement Parts



| 1. Reusable ventilation hose system with patient valveVentilation hose | WM 22520 WM 22647 | 5. Reusable PEEP valve infinitely adjustable from 2.5 to 10 mbar,22 mm external diameter | WM 3215 |
|--|----------------------------------|---|------------------------------|
| Patient valve, consisting of connection nozzle to patient hose, upper control element, bottom control element for spontaneous breathing and lip membrane | WM 3280 | 6. CapnoDura Combi disposable detector for endtidal CO2 measurerSet of 10 | WM 20760 ment WM 20770 |
| 2. Disposable ventilation hose system | WM 28110 | 7. Heat and moisture exchanging filter | WM 22162 |
| with patient valve | VVIVI ZOTTO | 8. Test bag | WM 1454 |
| Set of 10Set of 25Set of 50 | WM 15454 WM 15455 WM 15456 | 9. EasyLung for simulation of ventilation | WM 28625 |
| 3. Reusable hose protection sleeve | WM 8297 | 10. Battery | WM 28045 |
| 4. MEDUtrigger | WM 20900 | | |

for manual ventilation

Technical Data

| Dimensions incl. connections (WxHxD) | 100 x 145 x 90 mm | |
|--|---|--|
| Weight | 0.7 kg | |
| Temperature • For operation • For storage | -20 °C to +50 °C -40 °C to +70 °C | |
| Oxygen concentration | 100% | |
| Gas supply | Medical grade oxygen | |
| Oxygen consumption for device operation | 0 I/min | |
| Operating pressure | 2.7 bar to 6 bar | |
| Minimum input flow | 40 l/min | |
| Power supply | Lithium ion battery (3.6 V; 5.2 Ah) Battery runtime up to 100 h under standard condititions | |
| IPPV mode | | |
| Ventilation rate | Infinitely adjustable from 10 to 25 min ⁻¹ | |
| • Tidal volume (Vt) | Infinitely adjustable from 65 to 950 ml | |
| Correlation of tidal volume to ventilation frequency I:E ratio Max. ventilation pressure | f (min-1) 25 20 15 12 10 10 10 10 V, (ml) 65 100 150 300 500 600 800 950 1:2 or 1:3 20 or 45 mbar | |
| CPR mode | | |
| Tidal volume (Vt)Max. ventilation pressure | Adjustable from 65 ml to 950 ml 20 mbar or 45 mbar | |
| • I:E ratio | 1:1 | |
| Demand flow mode | Trigger 0.8 mbar Peak flow > = 40 l/min Switch-off pressure 3 mbar | |
| Visual and audio alarms | Airway pressure high Airway pressure low/Apnea Supply pressure low (2.7 bar O₂) Battery low | |
| Protection from dust and splashing water | IP54 | |
| Applied standards | DIN EN 794-3, ISO 10651-3, EN 1789, EN 60601-1, EN 60601-1-2, EN 60601-1-12, RTCA/DO-160 | |





All rights to design and specification modifications reserved.













Simply Professional

WEINMANN Emergency is a family-owned, internationally active medical technology company. With our mobile system solutions for emergency, transport and disaster medicine, we set standards for saving human lives. In close collaboration with professional users in emergency medical services, hospitals and military medical corps, we develop innovative medical products for ventilation and defibrillation. For more than 100 years we have offered our customers a high degree of reliability, extensive experience and quality made in Germany.

Headquarters

WEINMANN Emergency Medical Technology GmbH + Co. KG Frohbösestraße 12 22525 Hamburg Germany

T: +49 40 88 18 96-0

F: +49 40 88 18 96-480

T: +49 40 88 18 96-120 Customer Service

T: +49 40 88 18 96-122 After-Sales Service

E: info@weinmann-emt.de

Center for Production, Logistics, Service

WEINMANN Emergency Medical Technology GmbH + Co. KG Siebenstücken 14 24558 Henstedt-Ulzburg Germany

China

Weinmann (Shanghai) Medical Device Trading Co. Ltd. T: +86 21 52 30 22 25 ● info@weinmann-emt.cn

U.A.E.

WEINMANN Emergency Medical Technology GmbH + Co.KG (Branch) T: +971 432 100 31 ● info-dubai@weinmann-emt.com

France

WEINMANN Emergency France SARL – Paris – Les Ulis T: +33 1 69 41 51 20 ● info@weinmann-emt.fr

Russia

Weinmann SPb GmbH – St. Petersburg T: +7 812 633 30 82 • info@weinmann-emt.ru

Singapore

Weinmann Singapur PTE, Ltd.

T: +65 65 09 44 30 • info-singapore@weinmann-emt.sg

Spain

WEINMANN Emergency Medical Technology GmbH + Co. KG T: +34 91 79 01 137 • info-spain@weinmann-emt.es

USA

Weinmann Emergency LP

T: +1 770-274-2417 • info@weinmann-emergency.com





