



People. Health. Care.

Product Features at a Glance

Plaster of Paris (P.o.P)	Cellona® P.o.P.	Cellona® Xtra	Cellona® P.o.P.	Cellona® Xtra
	Bandage		Splint	
Fabric	17-thread count cotton gauze			
Product composition	α & β -calcium sulphate hemihydrate			
Number of Layer	roll	roll	folded in 4 layers	
Rigidity	very high			
Immersion Time	1 sec / meter			
Water Temperature	20 – 25 °C			
Setting Time	3 - 4 minutes	2.5 minutes	3 - 4 minutes	2.5 minutes
First Stability (safe for transport)	after 30 minutes			
Final Setting Time (allows full load bearing)	after 12 – 24 hours			
How Removed	Orthopaedic saw		Removable splint	
Used for	Circular cast		Splints & reinforce circular cast	

